



Copyright © 2003-2021 by Xoran Technologies, LLC. All Rights Reserved.

--- CONFIDENTIAL ---

License Server User Manual

Revision 1

Revision History

Rev	Date	Author	Project #	Description	DCO #
1	February 18, 2022	Slobodan Komazec	ECN 2019-030 / 5087	Initial release	2556

Contents

Revision History	2
Contents	3
Introduction.....	4
Log In	4
Create New License	5
Edit License.....	7
Associate License	7
Download License.....	8
Search License	8

License Server User Manual

Introduction

The purpose of this document is to show how this site can be used.

Log In

Enter username and password and press Log In. User must be member of “XORANTECH” domain.

Log In

Username

Password

Log In

© 2022 - Xoran Technologies, LLC.

The license form will be displayed.

Xoran License Server User Manual

Index

Create

Search: Clear

License ID	Computer ID	Institution	Application mode	Expiration Date	
5434534534		asfgasgsag	vTRON	2/15/2023	✕ 🔍 ⬆
4545		sassa	vTRON	2/11/2023	✕ 🔍 ⬆
NewLicense		New	MiniCAT	12/21/2022	✕ 🔍 ⬆
asdf	8CE3A2D2-E26B26A9-9E6D0D41-A3068068-E9C008A2	Xoran	MiniCAT	11/7/2022	✕ 🔍
asdfsadf1		asdfsadf1	MiniCAT	11/6/2022	✕ 🔍 ⬆
59	61	591	MiniCAT2D	11/4/2022	✕ 🔍
33333	62	ssss	MiniCAT	11/3/2022	✕ 🔍
99999	69	aaaa	MiniCAT2020	11/3/2022	✕ 🔍
76766767	70	Eee	VetCAT	11/3/2022	✕ 🔍
678	71	Set	MiniCAT	11/3/2022	✕ 🔍
12345	75	Lalaa	MiniCAT2020	11/3/2022	✕ 🔍

© 2022 - Xoran Technologies, LLC.

Create New License

Press the Create button in the upper left corner

Index



The screenshot shows a web interface with a header 'Index'. Below it, there is a button labeled 'Create' which is highlighted with a red rectangular box. Underneath the button, there are two labels: 'License ID' and 'Computer ID', each followed by a horizontal line representing an input field.

The new license entry form will be displayed. In General section choose whether it is a license, enter license id, institution, choose application mode and expiration date. Selecting the application mode may cause some options in the Details section to be checked.

General



The screenshot shows the 'General' section of a form. It contains the following fields and controls:

- Demo License**: A checkbox that is currently unchecked.
- License ID**: A text input field.
- Institution**: A text input field.
- Application Mode**: A dropdown menu with 'MiniCAT' selected.
- Expiration Date**: A text input field containing the value '2/16/2023 11:59:59 PM'.

In Details section choose License type, and, according to choosed type, some options could be checked. Press Create button to create a new license.

Details

License Type	<input type="text" value="v"/>
User Account Management	<input type="checkbox"/>
Patient Management	<input type="checkbox"/>
Acquisition	<input type="checkbox"/>
Reconstruction	<input type="checkbox"/>
3D	<input type="checkbox"/>
Modify study header	<input type="checkbox"/>
Delete Study	<input type="checkbox"/>
Study purge	<input type="checkbox"/>
Import Study	<input type="checkbox"/>
Export Study	<input type="checkbox"/>
Backup Push	<input type="checkbox"/>
Backup on Demand	<input type="checkbox"/>
Full Backup	<input type="checkbox"/>
Modify Actions	<input type="checkbox"/>
MWL Support	<input type="checkbox"/>
Radiographic Filters	<input type="checkbox"/>
Anisotropic Diffusion Filters	<input type="checkbox"/>
IQ Filters	<input type="checkbox"/>
MPR Position	<input type="checkbox"/>
Flip Horizontal	<input type="checkbox"/>

Create

[Back to List](#)

Edit License

To edit license, click on button for license edit, located on the right side of the screen:



A form similar to creating a license will appear. Enter the required changes and press the Update button.

Associate License

To associate license to specific computer id, click on button for license association, located on the right side of the screen:



Form for license association will appear:

Associate License

LicenseNumber

4545

ComputerID

Save

[Back to List](#)

© 2022 - Xoran Technologies, LLC.

Enter computer id and press Save button. After saving, button for license association will no longer be visible.

Download License

To download license, click on button for license download, located on the right side of the screen:



Search License

In order to search licenses whose License Id, Computer Id or institution contains a specific string, enter the string in the Search field and press Enter. To delete the search string, press the Clear button.

Search:

Institution	Application mode	Expiration Date
-------------	------------------	-----------------

To sort licenses by specific column, click on column name.